

# Detection of Cox2 in Formalin-Fixed, Paraffin-Embedded Rat and Mouse Tissue

## Reagents:

[1X Automation Buffer](#)  
[3% Hydrogen Peroxide](#)  
[Antibody Diluent](#)  
[Citrate Buffer](#)  
[DAB Chromagen](#)  
[Hematoxylin](#)

## Antibody Information:

Blocking Serum: Normal Goat Serum  
Jackson ImmunoResearch Laboratories, Inc.  
West Grove, PA 19390  
[www.jacksonimmuno.com](http://www.jacksonimmuno.com)  
1-800-367-5296  
Catalog #005-000-001

Avidin Biotin Blocking Kit  
Vector Laboratories, Inc.  
Burlingame, CA 94010  
[www.vectorlabs.com](http://www.vectorlabs.com)  
1-800-227-6666  
Catalog #SP-2001

Primary antibody: Rabbit anti-mouse Cox2  
Cayman Chemical  
Ann Arbor, MI 48108  
[www.caymanchem.com](http://www.caymanchem.com)  
1-800-364-9897  
Catalog #160106

Negative control serum: Normal Rabbit Serum  
Jackson ImmunoResearch Laboratories, Inc.  
West Grove, PA 19390  
[www.jacksonimmuno.com](http://www.jacksonimmuno.com)  
1-800-367-5296  
Concentration: 60 mg/ml  
Catalog #011-000-001

Secondary antibody: Biotinylated goat anti-rabbit IgG  
Vector Laboratories, Inc.  
Burlingame, CA 94010  
[www.vectorlabs.com](http://www.vectorlabs.com)  
1-800-227-6666  
Catalog #BA-1000

Label antibody: Vector EliteVectastain® ABC  
Vector Laboratories, Inc.  
Burlingame, CA 94010  
[www.vectorlabs.com](http://www.vectorlabs.com)  
1-800-227-6666  
Catalog #PK-6100

### **Staining Procedure:**

-Positive Control Tissue: Highest level of expression located at the distal vas deferens where it inserts into the bladder. Weak staining at the proximal end of the vas deferens.

-Stain localization: Peri-Nuclear/Cytoplasmic

Deparaffinize and hydrate slides through the following solutions.

Xylene	2 times	5 minutes
100% EtOH	2 times	3 minutes
95% EtOH	2 times	3 minutes
1X Automation Buffer	2 times	5 minutes

1. Quench endogenous peroxidase by placing slides in [3% hydrogen peroxide](#) for 15 minutes.
2. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.
3. Perform Heat Induced Epitope Retrieval using Microwave Oven.  
Place a full rack of slides in Tissue Tek™ container containing 200 mls 1X citrate buffer.  
MWO for 5 minutes at level 3  
Cool for 1 minute (Add 50 mls citrate buffer to container)  
MWO for 5 minutes at level 3. Temp after microwaving \_\_\_\_\_  
Cool 20 minutes at room temperature  
Rinse in distilled water 3 X 2 minutes each  
Place slides in buffer for 5 minutes
4. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

5. Block in 5% Normal Goat Serum for 20 minutes.

Lot#\_\_\_\_\_ Reconstituted Date\_\_\_\_\_

DO NOT RINSE SLIDES. CONTINUE TO AVIDIN-BIOTIN BLOCK.

6. Apply Avidin/Biotin block

Lot#\_\_\_\_\_ New Kit     yes   /   no

Apply avidin block - 15 min @ RT.

Quick rinse in 1X Automation Buffer.

Apply biotin block - 15 min @ RT.

No wash, wipe excess block and apply primary antibody

DO NOT RINSE SLIDES WITH BUFFER BEFORE ADDING PRIMARY ANTIBODY.

7. Apply primary antibody (Cox2) at 1:1000 dilution and incubate for one hour.

Lot#\_\_\_\_\_ Aliquoted   yes   /   no Date Aliquoted\_\_\_\_\_

Comment: The concentration of this antibody varies significantly from lot-to-lot. It is advisable to test the dilution (1:10, 1:100, 1:1000) of a new lot on a positive control before doing a large run.

For the negative control slides, match the protein concentration of the normal rabbit serum to the protein concentration of the primary antibody (Cox2) and use this to make the 1:1000 dilution. Apply to the slides and incubate for one hour.

Lot #\_\_\_\_\_ Reconstituted Date\_\_\_\_\_

8. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

9. Apply secondary antibody (Biotinylated goat anti-rabbit) at 1:500 dilution and incubate for 30 minutes.

Lot#\_\_\_\_\_ Reconstituted Date\_\_\_\_\_

10. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

11. Apply Label antibody and incubate for 30 minutes. (Prepare 30 minutes prior to use)

Exp. Date \_\_\_\_\_ New Kit     yes   /   no

12. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

13. Apply liquid Dako DAB Chromagen for 6 minutes in the dark.

(Add 1 drop of DAB per ml of substrate)

Lot #\_\_\_\_\_ Exp. Date\_\_\_\_\_ New Kit     yes   /   no

14. Rinse in tap water 3 minutes.

15. Counterstain with Modified Harris Hematoxylin for 30 seconds.
16. Rinse in tap water until water is clear.
17. Place slides in 1X Automation buffer for 1 minute with gentle agitation to blue slides.
18. Dehydrate through the following solutions.

95% Ethanol	1 change	3 minutes
100% EtOH	3 changes	3 minutes
Xylene	2 changes	5 minutes

19. Coverslip  
updated 01/14/2004