

Detection of CYP2C12 in Formalin-Fixed, Paraffin-Embedded Rat Tissue

Reagents:

[1X Automation Buffer](#)

[3% Hydrogen Peroxide](#)

[Antibody Diluent](#)

[Citrate Buffer](#)

[DAB Chromagen](#)

[Hematoxylin](#)

Antibody Information:

Blocking serum: Normal Goat Serum

Jackson ImmunoResearch Laboratories, Inc.

West Grove, PA 19390

www.jacksonimmuno.com

1-800-367-5296

Catalog # 005-000-121

Avidin Biotin Blocking Kit

Vector Laboratories, Inc.

Burlingame, CA 94010

www.vectorlabs.com

1-800-227-6666

Catalog #SP-2001

Primary antibody: Rabbit Polyclonal to Cytochrome P450 2C + 2C9 + 2C19 + 2C12

Abcam

Cambridge, MA 02139

www.abcam.com

1-888-772-2226

Catalog# ab22596-50

Negative control serum: Normal Rabbit Serum

Jackson ImmunoResearch Laboratories, Inc.

West Grove, PA 19390

www.jacksonimmuno.com

1-800-367-5296

Catalog #011-000-001

Secondary antibody: Biotinylated goat anti-rabbit IgG

Vector Laboratories, Inc.

Burlingame, CA 94010

www.vectorlabs.com

1-800-227-6666

Catalog # BA-1000

Label antibody: Peroxidase –conjugated Streptavidin SS Label

Biogenex

San Ramon, CA 94583

www.biogenex.com

1-800-421-4149

Catalog # HK330-9K

Staining Procedure

Positive Control Tissue: Rat Liver (upregulated)

Stain Localization: Centrilobular cytoplasmic staining pattern

Deparaffinize and hydrate slides through the following solutions:

Xylene	2 times	5 minutes
100% Ethanol	2 times	3 minutes
95% Ethanol	2 times	3 minutes
1X Automation Buffer	2 times	5 minutes

1. Quench endogenous peroxidase by placing slides in 3% hydrogen peroxide for 15 minutes.
2. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.
3. Unmasking Technique using the decloaker.
Add 500ml distilled water to the pan of the decloaker.
Place a full rack of slides in a Tissue Tek™ container containing 250ml of 1X citrate buffer solution.
Decloak for 5 minutes. Pressure _____
Depressurize for 10 minutes.
Remove pan top and cool for 10 minutes. Temperature after cooling _____
Rinse in distilled water two times for 3 minutes each.
4. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.
5. Block with 10% Normal Goat Serum for 20 minutes at room temperature.
Lot# _____ Reconstituted Date _____

DO NOT RINSE SLIDES. CONTINUE TO AVIDIN-BIOTIN.

6. Apply Avidin/Biotin block

Lot# _____ Exp. Date _____ New Kit: yes / no

Apply avidin block - 15 minutes at room temperature.

Quick rinse in 1X Automation Buffer

Apply biotin block - 15 minutes at room temperature.

Wipe excess block

DO NOT RINSE SLIDES WITH BUFFER BEFORE ADDING PRIMARY ANTIBODY.

7. Apply primary antibody (Cyp2C12) at 1:1000 dilution and incubate for one hour at room temperature.

Lot# _____ Date Aliquoted _____

For negative control slides, normalize the normal rabbit serum to match the protein concentration of the primary antibody (Cyp2C12), and use this to make a 1:1000 dilution. Apply to the slides and incubate for one hour at room temperature.

Lot# _____ Reconstituted Date _____

8. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

9. Apply Goat anti-rabbit secondary antibody at a 1:500 dilution and incubate for 30 minutes at room temperature.

Lot# _____ Reconstituted Date _____

10. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

11. Apply the Biogenex Supersensitive Label and incubate 30 minutes at room temperature.

Lot# _____ Exp. Date _____

12. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

13. Apply liquid Dako DAB Chromagen for 6 minutes in the dark.

(Add 1 drop of DAB per ml of substrate)

Lot# _____ Exp. Date _____ New Kit: yes / no

14. Rinse in tap water 3 minutes.

15. Counterstain with Modified Harris Hematoxylin for 20 seconds.

16. Rinse in tap water until water is clear.

17. Gently agitate slides in 1X Automation Buffer until blue.

18. Dehydrate through the following solutions.

95% Ethanol	1 change	3 minutes
100% Ethanol	3 changes	3 minutes
Xylene	2 changes	5 minutes

19. Coverslip

Updated 08/21/06