

Advances in Microarray Hybridizations at NIEHS

Sherry Grissom

Microarray User's Meeting

June 21, 2002

Evolution of the Microarray Lab

Where We Were	Where We Are	Why It's Good
200 μ g RNA needed per microarray chip	100-200 μ g RNA/ microarray experiment (4 chips)	DIR likes us better!
Microcon cleanup	Qiagen cleanup	Time, clean chips
96-well PCRs by hand	Qiagen BioRobot	Saves time
Single Manual Sequencing	Qiagen BioRobot (96-well)	Saves time
4 & 8 pin printing	32 pin printing	Saves time
Axon Scanner	Agilent Scanner	Time to do other things, not subjective
ArraySuite 1.3	ArraySuite 2.0	Quality & Signal-to-Noise Cutoffs

Indirect vs. Direct Labeling

Protocol	µg RNA	Method	Time
NIEHS	25–50	Direct incorporation of dTTP analogs Cy3-dUTP & Cy5-dUTP at RT step	5-6 hours
3DNA (Genisphere)	1	Indirect, Primer oligo with capture sequence attaches at dT tail, bind capture reagent that recognizes capture sequence	5-6 hours
Atlas Glass (Clontech)	10-20	Indirect RT incorporation of aminoallyl-dUTP, Amine reactive dye added.	8 hours

Cy3 mean=1881
Background=265

Cy5 mean=1551
Background=441





Image Name: p54s75_635_flpmrk.ip
p54s75_532_flp.ip

☐ Histogram ☒ Scatter Plot



Data from:

- ☒ All targets
☐ Control targets

Calibration by:

- ☒ Internal Controls
☐ All targets
☐ Background

Calibration Method:

- ☒ Ratio Distribution
☐ Log-Normal

☒ Log Scale

☒ Calibrated Result

Confidence Level: 99.00

Intensity From: 200

Ratio Lower Limit: 0.76

Intensity To: 65535

Ratio Upper Limit: 1.32

Target Size From: 100

Refine Stats

CV = 0.075

M = 0.901

Default Setting

Ratio Outliers

Update dataSheet

Equal Intensities

Exit

Refresh



TargetLocator, 1.3, © 1997, NIH/NHGRI/LCG

- Very tight scatterplot with few outliers
- Low variation and M=1 reflect a balance between channels.

Direct Labeling by NIEHS SOP

Cy3 mean=303
Background=101

Cy5 mean=93
Background=80



- Poor signal-to-noise ratio.
- Wide scatterplot reflects poor dye incorporation.
- High CV and M value reflect imbalance.

3DNA by Genisphere

Product Evaluations

- We have also looked at Molecular Probes ARES, Stratagene Fair Play, other reverse transcriptases and Agilent Human 2 cDNA microarrays
- At this time, our current direct labeling protocol is best for in-house arrays.
- Genisphere 3DNA kit is constantly being updated, will try again.
- Soon will retry Clontech's Atlas Array protocol (been tried before with promising results) and Invitrogen's Fluoroscript labeling kit.

**How the Agilent Bioanalyzer
has made my life a bit easier
and how it can help yours**

Utilities of the Agilent Bioanalyzer



RNA Quality Measure: qualitative tool to measure RNA integrity for labeling.

- Degradation of either 18S or 28S rRNA
- Concentration
- Sometimes, genomic DNA contamination

Bioanalyzer Test Submissions

Bioanalyzer Test Submission Form

This form must be included with all samples submitted to the NIEHS Microarray Center for test on the Agilent Bioanalyzer.

Name and Contact Information (room #, phone #):

Email Address:

Principal Investigator:

Project (proposal) Title:

Sample Submission Requirements

(samples will be returned if all criteria are not met):

Please check boxes when criteria are met

- ☐ a. RNA concentration of a 100-200 ng/ul RNA in nuclease-free water.
- ☐ b. Minimum volume of 5 ul in an autochved tube, delivered on dry ice.
- ☐ c. If two or more samples are isolated with different procedures/protocols please indicate the differences (this could help us give better advice in the future).
- ☐ d. Did someone in the Microarray Center give any advice/help on your isolation, if so, please indicate this and the member of the Center who did.

*Please provide unambiguous names for RNA (no 3-letter codes or numbering 1-12). If we find that the annotations given to separate samples are confusing, we will return this form for further clarification.

**Please try to limit test samples to no more than 12 at a time.

***Because they are test samples, RNA samples submitted in this way are not considered a priority for the Bioanalyzer. They will be run as soon as space is available.

For Bioanalyzer test sample submission, please contact Sherry Grissom either by phone (x0747) or email (grissom2).

RNA INFORMATION

Tube #	RNA Name	RNA Type (tissue, cell line)	RNA Species (total, mRNA)	Prep Date	Prep Procedure
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					

<http://dir.niehs.nih.gov/microarray/methods/bioanalyzer.pdf>

Utilities of the Agilent Bioanalyzer




Labeled RNA Quality Measure: qualitative tool to measure labeled RNA integrity for hybridization.

- Is there a problem with the labeling and reverse transcription reaction?
- “Cheaper” way to optimize protocol
- Others use spectrophotometer and mini-gels to test the integrity of their labeled probe.


Spectrophotometer vs. Agilent Bioanalyzer

- Spectrophotometer output is absorbance at a given wavelength.
- From the absorbance measurements, can calculate concentration, labeling density, and frequency of dye incorporation.
- No way to tell about dye carry over.
- Need at least 50 μL sample.
- Bioanalyzer output is in form of an electropherogram.
- Directly outputs **semi-quantitative** labeled RNA concentrations and unincorporated Cy5 amounts.
- Gives relative idea of amount of dye being carried over.
- Only need 1 μL of sample


Dilution Scheme





Comparison of Spectrophotometer and Bioanalyzer



Cy5 Incorporation as a Measure of RT reaction efficiency



NMC's QC RNA
(we know it labels)



Questionable RNA
(we know it doesn't label)



Summer Exp. 1- Negative Control Background Subtraction

- Should we be doing global or local background subtraction? Currently, we are assuming that what is “beside” a spot is equal to what is “underneath” that spot. Is this a correct assumption?
- Print a rat chip, replace 566 spots on the array (17 per sub-array) with negative control cDNA spots (Arabidopsis Cab), perform hybridizations. Data analysis performed by Terry Speed at UC-Berkeley.
- Theoretically, there should be no signal obtained from this spot because there will be no cross-hybridization. In practice, signal will be obtained from these spots as a result of cross-hybridization and autofluorescence.

Summer Exp. 2 -Qualitymetrics

- Artificially derived sequences from yeast intergenic regions (YIRs) with no cross-reactivity among species.
- Spike-In = YIR cDNA spotted and corresponding mRNA sequence of known concentration “spiked” into reverse transcription reaction.
 - Calibration controls are used to determine the dynamic range of the reaction.
 - Ratio controls tell us how well each of the labeling reactions went.
- Negative Controls
 - No cross-hybridization (should be the same as no spot)
 - How do they compare to buffer blanks?
 - Is there uniform distribution across the chip?
- Positive Controls
 - Was the label a success?
 - Is there uniform distribution across the chip?

A few intriguing questions for later...

- How low can you go? Getting the maximum amount of signal for the minimum quantities of RNA and other reagents.
- How can we measure genomic contamination, why is it detrimental to the microarray process and what can we do about it?

Future Goals

- Bioanalyzer for measuring Cy5 incorporation to:
 - Can we use less RNA and get comparable signals?
 - Optimize RT enzyme for time and temperature.
 - Optimize Cy dye concentrations to reflect lower amounts of RNA used for hybridizations.
- Continue DIR projects (since August 2001, I have performed ~400 hybridizations).

Thanks

- Jeff Tucker
- Cindy Afshari
- Rick Paules
- Astrid Haugen
- Danica Ducharme
- Valerie Moorman
- Neysa Garner



For questions, please visit <http://dir.niehs.nih.gov/microarray>
or D226